

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 210

Grad. pipettes, 2:0,1 ml, AS, conf. cert.

Batch: 04-21




Serial:

Measuring Nr. 44138

Batch size:

225

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	1,9974 ml	1,2004 ml	0,4012 ml
Deviation:	0,0028 ml	0,0021 ml	0,0016 ml
Precision (CV):	0,1402 %	0,1780 %	0,3935 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

03.02.2021

Balance:

Sartorius -
LE244S-OCE

DE-M