

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 515 510 BL

**Grad. pipette, 5 : 0,1 ml, blue grad., AS, conf. cert.**




**Batch:** 19-19

**Serial:**

**Measuring Nr.** 29882

**Batch size:** 150

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	5,0029 ml	3,0033 ml	1,0057 ml
Deviation:	0,0048 ml	0,0045 ml	0,0028 ml
Precision (CV):	0,0962 %	0,1492 %	0,2841 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	14.05.2019
Balance:	Sartorius - LE244S-OCE

DE-M