

Certificate of performance



Article: 3 515 510

Grad. pipettes, 5:0,1 ml, AS, conf. cert.




Batch: 21-18

Serial:

Measuring Nr.: 21423

Batch size: 280

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	4,9934 ml	2,9930 ml	0,9978 ml
Deviation:	0,0080 ml	0,0048 ml	0,0021 ml
Precision (CV):	0,1597 %	0,1616 %	0,2140 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	44
Date:	04.06.2018
Balance:	Sartorius - LE244S-OCE

