

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 710

Grad. pipettes, 10:0,1 ml, AS, conf. cert.




Batch: 05-23

Serial:

Measuring Nr. 61671

Batch size: 500

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	10,0127 ml	6,0093 ml	2,0062 ml
Deviation:	0,0103 ml	0,0081 ml	0,0040 ml
Precision (CV):	0,1027 %	0,1343 %	0,2007 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	13.02.2023
Balance:	Sartorius - LE244S-OCE

DE-M