

# Certificate of performance



**Article:** 3 515 810

**Grad. pipettes, 20:0,1 ml, AS, conf. cert.**




**Batch:** 36-18

**Serial:**

**Measuring Nr.** 23739

**Batch size:** 280

Volume:	20,000 ml	12,000 ml	4,000 ml
<b>Tolerance:</b>	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
<b>Average value:</b>	19,9615 ml	11,9671 ml	3,9892 ml
<b>Deviation:</b>	0,0217 ml	0,0129 ml	0,0094 ml
<b>Precision (CV):</b>	0,1086 %	0,1074 %	0,2350 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 9

**Date:** 04.09.2018

**Balance:** Sartorius - LE244S-OCE

