

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 515 810

**Grad. pipettes, 20:0,1 ml, AS, conf. cert.**




**Batch:** 34-23

**Serial:**

**Measuring Nr.** 64372

**Batch size:** 150

Volume:	20,000 ml	12,000 ml	4,000 ml
<b>Tolerance:</b>	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
<b>Average value:</b>	19,9963 ml	11,9963 ml	4,0031 ml
<b>Deviation:</b>	0,0246 ml	0,0125 ml	0,0076 ml
<b>Precision (CV):</b>	0,1231 %	0,1045 %	0,1901 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	07.09.2023
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

DE-M