

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 515 910

**Graduated pipettes class AS 25 ml**




**Batch:** 45-20

**Serial:**

**Measuring Nr.:** 42331

**Batch size:** 1200

Volume:	25,000 ml	15,000 ml	5,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	25,0066 ml	15,0027 ml	5,0004 ml
Deviation:	0,0246 ml	0,0198 ml	0,0111 ml
Precision (CV):	0,0985 %	0,1321 %	0,2212 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	10.11.2020
Balance:	Sartorius - LE244S-OCE

