

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 011

Enzyme test pipette, 0,1: 0,001, AS, conf. cert.




Batch: 06-19

Serial:

Measuring Nr. 27510

Batch size: 280

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	± 0,0030 ml	± 0,0030 ml	± 0,0030 ml
Average value:	0,0990 ml	0,0586 ml	0,0185 ml
Deviation:	0,0006 ml	0,0007 ml	0,0006 ml
Precision (CV):	0,5655 %	1,1123 %	3,0875 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	05.02.2019
Balance:	Sartorius - LE244S-OCE

