

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 517 011

**Enzyme test-pipette, 0,1:0,001, AS, conf. cert.**

**Batch:** 09-20




**Serial:**

**Measuring Nr.** 36506

**Batch size:**

305

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	± 0,0030 ml	± 0,0030 ml	± 0,0030 ml
Average value:	0,1000 ml	0,0598 ml	0,0193 ml
Deviation:	0,0008 ml	0,0005 ml	0,0009 ml
Precision (CV):	0,8184 %	0,8249 %	4,2763 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	05.03.2020
Balance:	Sartorius - MSE224S-0CE -DU

