

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 011

Pipette f. encym. analysis., 0,1:0,001 ml, AS, conf. cert.




Batch: 27-21

Serial:

Measuring Nr.: 47794

Batch size: 482

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	± 0,0030 ml	± 0,0030 ml	± 0,0030 ml
Average value:	0,0996 ml	0,0588 ml	0,0188 ml
Deviation:	0,0007 ml	0,0007 ml	0,0007 ml
Precision (CV):	0,7006 %	1,2259 %	3,6593 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	12.07.2021
Balance:	Sartorius - LE244S-OCE

