

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 517 021

**Pipette f. encym. analysis., 0,2:0,001 ml, AS, conf. cert.**




**Batch:** 35-21

**Serial:**

**Measuring Nr.** 49345

**Batch size:** 150

Volume:	0,2000 ml	0,1200 ml	0,0400 ml
<b>Tolerance:</b>	± 0,0030 ml	± 0,0030 ml	± 0,0030 ml
<b>Average value:</b>	0,2001 ml	0,1197 ml	0,0390 ml
<b>Deviation:</b>	0,0008 ml	0,0004 ml	0,0007 ml
<b>Precision (CV):</b>	0,4067 %	0,3518 %	1,6373 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	28.09.2021
<b>Balance:</b>	Sartorius - LE244S-OCE

DE-M