

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 517 051 BL

**Pipette f.encym.analysis., 0,5:0,01 ml, AS, conf. cert.**




**Batch:** 25-19

**Serial:**

**Measuring Nr.** 31194

**Batch size:** 36

Volume:	0,5000 ml	0,3000 ml	0,1000 ml
<b>Tolerance:</b>	± 0,0050 ml	± 0,0050 ml	± 0,0050 ml
<b>Average value:</b>	0,4989 ml	0,2999 ml	0,0986 ml
<b>Deviation:</b>	0,0009 ml	0,0012 ml	0,0008 ml
<b>Precision (CV):</b>	0,1770 %	0,4138 %	0,8215 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12 699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	03.07.2019
<b>Balance:</b>	Sartorius - LE244S-OCE

