

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 051

Pipettes f. encym. analysis, 0,5:0,01 ml, AS, conf. cert.




Batch: 32-23

Serial:

Measuring Nr. 63967

Batch size: 280

Volume:	0,5000 ml	0,3000 ml	0,1000 ml
Tolerance:	± 0,0050 ml	± 0,0050 ml	± 0,0050 ml
Average value:	0,4983 ml	0,2985 ml	0,0992 ml
Deviation:	0,0010 ml	0,0010 ml	0,0006 ml
Precision (CV):	0,1957 %	0,3188 %	0,5565 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12 699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	15.08.2023
Balance:	Sartorius - LE244S-OCE

