

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 101 BL

Pipettes f. enzym. analysis, 1:0,01 ml, AS, conf. cert.




Batch: 04-18

Serial: -

Measuring Nr.: 17756

Batch size: 72

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	1,0006 ml	0,5994 ml	0,1994 ml
Deviation:	0,0021 ml	0,0015 ml	0,0005 ml
Precision (CV):	0,2063 %	0,2522 %	0,2581 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	01.02.2018
Balance:	Sartorius - LE244S-OCE

