

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 101

Pipettes f. encym. analysis, 1:0,01 ml, AS, conf. cert.




Batch: 14-20

Serial:

Measuring Nr. 37419

Batch size: 143

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	1,0016 ml	0,5996 ml	0,1997 ml
Deviation:	0,0016 ml	0,0011 ml	0,0007 ml
Precision (CV):	0,1607 %	0,1830 %	0,3726 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	24.04.2020
Balance:	Sartorius - LE244S-OCE

