

# Certificate of performance



**Article:** 3 517 101

**Pipettes f. encym. analysis, 1:0,01 ml, AS, conf. cert.**




**Batch:** 35-23

**Serial:**

**Measuring Nr.** 64412

**Batch size:** 185

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	0,9999 ml	0,6003 ml	0,2007 ml
Deviation:	0,0024 ml	0,0015 ml	0,0012 ml
Precision (CV):	0,2370 %	0,2575 %	0,5873 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	14.09.2023
Balance:	Sartorius - LE244S-OCE

