

# Certificate of performance



**Article:** 3 517 202

**Pipettes f. encym. analysis, 2:0,02 ml, AS, conf. cert.**




**Batch:** 24-22

**Serial:**

**Measuring Nr.** 56579

**Batch size:** 128

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	2,0037 ml	1,2027 ml	0,4042 ml
Deviation:	0,0022 ml	0,0018 ml	0,0009 ml
Precision (CV):	0,1087 %	0,1512 %	0,2149 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	20.06.2022
Balance:	Sartorius - LE244S-OCE

