

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 202

Pipettes f. encym. analysis, 2:0,02 ml, AS, conf. cert.




Batch: 10-24

Serial:

Measuring Nr. 66313

Batch size: 105

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	2,0018 ml	1,2033 ml	0,4022 ml
Deviation:	0,0029 ml	0,0031 ml	0,0015 ml
Precision (CV):	0,1453 %	0,2555 %	0,3634 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	12.03.2024
Balance:	Sartorius - LE244S-OCE

DE-M