

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 302

Pipettes f. encym. analysis, 3,0:0,02 ml, AS, conf. cert.




Batch: 17-18

Serial: -

Measuring Nr.: 20869

Batch size: 74

Volume:	3,0000 ml	1,8000 ml	0,6000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	2,9997 ml	1,8015 ml	0,6019 ml
Deviation:	0,0042 ml	0,0029 ml	0,0010 ml
Precision (CV):	0,1400 %	0,1631 %	0,1621 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	44
Date:	22.05.2018
Balance:	Sartorius - LE244S-OCE

