

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 505

Pipettes f. encym. analysis, 5:0,05 ml, AS, conf. cert.




Batch: 26-18

Serial:

Measuring Nr.: 21864

Batch size: 280

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	4,9979 ml	2,9981 ml	1,0016 ml
Deviation:	0,0055 ml	0,0030 ml	0,0015 ml
Precision (CV):	0,1099 %	0,1014 %	0,1467 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	44
Date:	02.07.2018
Balance:	Sartorius - LE244S-OCE

