

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 517 505

**Pipettes for enzym. analysis, 5:0,05 ml, AS, conf. cert.**




**Batch:** 40-22

**Serial:** -

**Measuring Nr.:** 59781

**Batch size:** 50

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
<b>Tolerance:</b>	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
<b>Average value:</b>	4,9972 ml	2,9973 ml	0,9995 ml
<b>Deviation:</b>	0,0061 ml	0,0048 ml	0,0016 ml
<b>Precision (CV):</b>	0,1228 %	0,1605 %	0,1600 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

24.10.2022

**Balance:**

Sartorius -  
MSE224S-0CE  
-DU

DE-M