

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 710

Grad. pipettes for enz. analysis, 10 ml, AS, conf. cert.

Batch: 26-17




Serial: -

Measuring Nr.: 11731

Batch size:

280

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	10,0068 ml	6,0004 ml	1,9998 ml
Deviation:	0,0099 ml	0,0053 ml	0,0028 ml
Precision (CV):	0,0988 %	0,0875 %	0,1386 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

03.07.2017

Balance:

Sartorius -
LE244S-OCE

DE-M