

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 710

Grad. pipettes for enz. analysis, 10 ml, AS, conf. cert.




Batch: 08-18

Serial:

Measuring Nr. 19188

Batch size: 239

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9889 ml	5,9971 ml	2,0026 ml
Deviation:	0,0113 ml	0,0063 ml	0,0032 ml
Precision (CV):	0,1126 %	0,1046 %	0,1608 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	44
Date:	27.03.2018
Balance:	Sartorius - LE244S-OCE

