

Certificate of performance



Article: 3 517 710

Pipettes f. encym. analysis, 10:0,1 ml, AS, conf. cert.




Batch: 49-19

Serial:

Measuring Nr. 35100

Batch size: 280

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9951 ml	5,9973 ml	2,0069 ml
Deviation:	0,0109 ml	0,0076 ml	0,0048 ml
Precision (CV):	0,1092 %	0,1275 %	0,2403 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	11.12.2019
Balance:	Sartorius - LE244S-OCE

