

# Certificate of performance



**Article:** 3 517 710

**Pipettes f. encym. analysis, 10:0,1 ml, AS, conf. cert.**




**Batch:** 17-21

**Serial:**

**Measuring Nr.** 46481

**Batch size:** 700

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9957 ml	5,9995 ml	2,0038 ml
Deviation:	0,0152 ml	0,0112 ml	0,0057 ml
Precision (CV):	0,1516 %	0,1866 %	0,2839 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	05.05.2021
Balance:	Sartorius - LE244S-OCE

