

Certificate of performance

Article: 3 518 011 BL

Grad. pipettes, 0,1 ml, Schellbach, AS, conf. cert.




Batch: 12-20

Serial:

Measuring Nr. 42835

Batch size: 5

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	0,1013 ml	0,0608 ml	0,0196 ml
Deviation:	0,0018 ml	0,0003 ml	0,0003 ml
Precision (CV):	1,8315 %	0,4821 %	1,2633 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	77
Date:	25.11.2020
Balance:	Sartorius - MSE224S-OCE -DU