

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 518 101 BL

**Grad. pipettes, 1,0:0,01 ml, Schellbach, AS, conf. cert.**




**Batch:** 49-20

**Serial:**

**Measuring Nr.** 43239

**Batch size:** 147

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
<b>Tolerance:</b>	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
<b>Average value:</b>	0,9994 ml	0,5989 ml	0,1990 ml
<b>Deviation:</b>	0,0017 ml	0,0017 ml	0,0016 ml
<b>Precision (CV):</b>	0,1706 %	0,2847 %	0,8187 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	12.12.2020
<b>Balance:</b>	Sartorius - LE244S-OCE

