

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 518 710 BL

Grad. pipette, 10 ml, Schellbach, AS, conf. cert.




Batch: 18-19

Serial:

Measuring Nr.: 30125

Batch size: 150

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9959 ml	5,9961 ml	1,9977 ml
Deviation:	0,0130 ml	0,0053 ml	0,0076 ml
Precision (CV):	0,1299 %	0,0884 %	0,3806 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	23.05.2019
Balance:	Sartorius - LE244S-OCE

