

Certificate of performance



Article: 3 518 710 BL

Grad. pipette, 10 ml, Schellbach, AS, conf. cert.




Batch: 08-21

Serial:

Measuring Nr. 44718

Batch size: 150

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9677 ml	5,9873 ml	1,9942 ml
Deviation:	0,0057 ml	0,0101 ml	0,0054 ml
Precision (CV):	0,0568 %	0,1688 %	0,2692 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	02.03.2021
Balance:	Sartorius - LE244S-OCE

