

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 518 710 BL

Grad. pipette, 10 ml, Schellbach, AS, conf. cert.




Batch: 10-24

Serial:

Measuring Nr. 66444

Batch size: 150

| Volume: | 10,000 ml | 6,0000 ml | 2,0000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0500 ml | ± 0,0500 ml | ± 0,0500 ml |
| Average value: | 10,0076 ml | 5,9865 ml | 1,9955 ml |
| Deviation: | 0,0105 ml | 0,0077 ml | 0,0065 ml |
| Precision (CV): | 0,1051 % | 0,1280 % | 0,3268 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 19.03.2024 |
| Balance: | Sartorius - LE244S-OCE |

DE-M