

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 040 BL

Vol. pipettes, 4 ml, blue grad., AS, conf. cert.

Batch: 47-24


Serial:

Measuring Nr. 71307

Batch size:

257

| | |
|------------------------|-----------------|
| Volume: | 4,0000 ml |
| Tolerance: | $\pm 0,0150$ ml |
| Average value: | 4,0024 ml |
| Deviation: | 0,0037 ml |
| Precision (CV): | 0,0931 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

13.01.2025

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M