

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 040 BL USP

Vol. pipettes, 4 ml, blue grad., USP, AS, conf. cert.,

Batch: 16-24


Serial:

Measuring Nr. 66890

Batch size:

60

Volume:	4,0000 ml
Tolerance:	$\pm 0,0100$ ml
Average value:	4,0019 ml
Deviation:	0,0015 ml
Precision (CV):	0,0384 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

22.04.2024

Balance:

Sartorius -
LE244S-OCE

DE-M