

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 535 040 BL

**Vol. pipettes, 4 ml, blue grad., AS, conf. cert.**


**Batch:** 11-19

**Serial:**

**Measuring Nr.** 28956

**Batch size:** 500

|                        |                  |
|------------------------|------------------|
| <b>Volume:</b>         | <b>4,0000 ml</b> |
| <b>Tolerance:</b>      | $\pm 0,0150$ ml  |
| <b>Average value:</b>  | 4,0045 ml        |
| <b>Deviation:</b>      | 0,0023 ml        |
| <b>Precision (CV):</b> | 0,0582 %         |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                           |
|------------------|---------------------------|
| <b>Operator:</b> | 0                         |
| <b>Date:</b>     | 02.04.2019                |
| <b>Balance:</b>  | Sartorius -<br>LE244S-OCE |

DE-M