

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 0585

Volumetric pipettes, 5,85 ml, AS, conf. cert

Batch: 24-20


Serial:

Measuring Nr. 38999

Batch size:

60

| | |
|------------------------|------------------|
| Volume: | 5,8500 ml |
| Tolerance: | $\pm 0,0150$ ml |
| Average value: | 5,8455 ml |
| Deviation: | 0,0020 ml |
| Precision (CV): | 0,0347 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

06.07.2020

Balance:

Sartorius -
LE244S-OCE

DE-M