

Certificate of performance

Article: 3 535 060 BL

Volumetric pipettes, 6 ml, blue grad., AS, conf. cert.


Batch: 22-20

Serial:

Measuring Nr. 38709

Batch size: 150

Volume:	6,0000 ml
Tolerance:	$\pm 0,0150$ ml
Average value:	5,9978 ml
Deviation:	0,0039 ml
Precision (CV):	0,0642 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

16.06.2020

Balance:

Sartorius -
LE244S-0CE