

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 080 USP

Volumetric pipettes, 8 ml, USP, AS, conf. cert.


Batch: 03-21

Serial:

Measuring Nr. 44616

Batch size: 90

Volume:	8,0000 ml
Tolerance:	± 0,0200 ml
Average value:	7,9996 ml
Deviation:	0,0033 ml
Precision (CV):	0,0415 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	22.02.2021
Balance:	Sartorius - LE244S-OCE

DE-M