

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 200 BL

Volumetric pipettes, 20 ml, blue grad., AS, conf. cert.

Batch: 01-21


Serial:

Measuring Nr. 44598

Batch size:

279

| | |
|------------------------|-----------------|
| Volume: | 20,000 ml |
| Tolerance: | $\pm 0,0300$ ml |
| Average value: | 20,0033 ml |
| Deviation: | 0,0039 ml |
| Precision (CV): | 0,0196 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 18.02.2021 |
| Balance: | Sartorius - LE244S-OCE |

DE-M