

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 320

Volumetric pipette, 32 ml, AS, conf. cert.


Batch: 27-19

Serial:

Measuring Nr. 33615

Batch size: 13

Volume:	32,000 ml
Tolerance:	$\pm 0,0300$ ml
Average value:	32,0211 ml
Deviation:	0,0025 ml
Precision (CV):	0,0078 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	11.10.2019
Balance:	Sartorius - MSE224S-0CE -DU

DE-M