

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 400 BL USP

Volumetric pipettes, 40 ml, blue, USP, AS, conf.cert,


Batch: 03-21

Serial:

Measuring Nr. 44576

Batch size: 50

Volume:	40,0000 ml
Tolerance:	$\pm 0,0500$ ml
Average value:	39,9827 ml
Deviation:	0,0051 ml
Precision (CV):	0,0126 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	18.02.2021
Balance:	Sartorius - LE244S-OCE

DE-M