

# Certificate of performance



ISO 9001  
CERTIFIED

Article: 3 535 400 BL

Volumetric pipettes, 40 ml, blue grad., AS, conf. cert.

Batch: 04-18


Serial:

Measuring Nr. 17755

Batch size:

180

Volume:	40,0000 ml
Tolerance:	$\pm 0,0500$ ml
Average value:	39,9991 ml
Deviation:	0,0088 ml
Precision (CV):	0,0220 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	31.01.2018
Balance:	Sartorius - LE244S-0CE

DE-M