

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 535 500 USP

**Vol. pipette, 50 ml, USP, AS, conf. cert**


**Batch:** 29-23

**Serial:**

**Measuring Nr.** 63890

**Batch size:** 90

<b>Volume:</b>	50,000 ml
<b>Tolerance:</b>	± 0,0500 ml
<b>Average value:</b>	49,9957 ml
<b>Deviation:</b>	0,0131 ml
<b>Precision (CV):</b>	0,0262 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 7

**Date:** 08.08.2023

**Balance:** Sartorius -  
LE244S-OCE

DE-M