

Certificate of performance

Article: 3 535 900 BL

Volumetric pipettes, 100 ml, blue grad., AS, conf. cert.

Batch: 49-19


Serial:

Measuring Nr. 35095

Batch size:

500

| | |
|------------------------|-------------------|
| Volume: | 100,000 ml |
| Tolerance: | $\pm 0,0800$ ml |
| Average value: | 100,0205 ml |
| Deviation: | 0,0209 ml |
| Precision (CV): | 0,0209 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

09.12.2019

Balance:

Sartorius -
MSE224S-0CE
-DU