

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 103 B

Burettes, 10 ml, str. stopc., amber grad., AS, conf. cert.




Batch: 16-23

Serial:

Measuring Nr. 63151

Batch size: 95

| Volume: | 10,000 ml | 6,0000 ml | 2,0000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0200 ml | ± 0,0200 ml | ± 0,0200 ml |
| Average value: | 9,9962 ml | 5,9978 ml | 2,0054 ml |
| Deviation: | 0,0064 ml | 0,0074 ml | 0,0039 ml |
| Precision (CV): | 0,0636 % | 0,1227 % | 0,1934 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 26.05.2023 |
| Balance: | Sartorius - LE244S-OCE |

DE-M