

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 103

Burettes, 10 ml, straight stopc., PTFE, AS, conf. cert.




Batch: 05-23

Serial:

Measuring Nr. 62915

Batch size: 50

Volume:	10,000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	9,9944 ml	5,9935 ml	2,0004 ml
Deviation:	0,0051 ml	0,0058 ml	0,0013 ml
Precision (CV):	0,0510 %	0,0970 %	0,0646 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 7

Date: 09.05.2023

Balance: Sartorius -
LE244S-OCE

DE-M