

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 010 502

**Burettes, 50 ml, str. stopc., solid plug, AS, conf. cert.**




**Batch:** 46-20

**Serial:**

**Measuring Nr.** 70016

**Batch size:** 5

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0027 ml	29,9723 ml	9,9885 ml
Deviation:	0,0155 ml	0,0039 ml	0,0107 ml
Precision (CV):	0,0310 %	0,0131 %	0,1067 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	77
Date:	05.11.2024
Balance:	Sartorius - LE244S-OCE

