

Certificate of performance



Article: 3 010 503 B

Burettes, 50 ml, str. stopc., amber grad., AS, conf. cert.




Batch: 16-24

Serial:

Measuring Nr. 67693

Batch size: 600

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9987 ml	29,9865 ml	10,0008 ml
Deviation:	0,0137 ml	0,0108 ml	0,0094 ml
Precision (CV):	0,0274 %	0,0361 %	0,0945 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	14.06.2024
Balance:	Sartorius - LE244S-OCE

