

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 504

Burettes 50 ml, lat. stopc., glass-plug, AS, conf. cert.




Batch: 46-20

Serial:

Measuring Nr. 64398

Batch size: 5

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0035 ml	29,9764 ml	9,9899 ml
Deviation:	0,0016 ml	0,0062 ml	0,0011 ml
Precision (CV):	0,0031 %	0,0206 %	0,0111 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	14.09.2023
Balance:	Sartorius - LE244S-OCE

DE-M