

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 053

Burettes, 5 ml, str. stopc., PTFE- plug, AS, conf. cert.




Batch: 40-22

Serial:

Measuring Nr. 61421

Batch size: 49

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	4,9958 ml	2,9965 ml	1,0031 ml
Deviation:	0,0022 ml	0,0027 ml	0,0023 ml
Precision (CV):	0,0436 %	0,0904 %	0,2282 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	31.01.2023
Balance:	Sartorius - LE244S-OCE

