

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 253

**Burettes, 25 ml, str. stopc., PTFE plug, AS, conf. cert.**




**Batch:** 22-22

**Serial:**

**Measuring Nr.** 61928

**Batch size:** 49

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9993 ml	14,9970 ml	5,0057 ml
Deviation:	0,0059 ml	0,0097 ml	0,0060 ml
Precision (CV):	0,0237 %	0,0644 %	0,1191 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

14.03.2023

**Balance:**

Sartorius -  
LE244S-OCE

DE-M