

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 503

**Burettes, 50 ml, str. stopc., PTFE plug, AS, conf. cert.**




**Batch:** 10-24

**Serial:**

**Measuring Nr.** 69450

**Batch size:** 72

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9903 ml	29,9823 ml	9,9843 ml
Deviation:	0,0169 ml	0,0103 ml	0,0119 ml
Precision (CV):	0,0339 %	0,0344 %	0,1189 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	09.10.2024
Balance:	Sartorius - MSE224S-0CE -DU

